


Product Datasheet

Human C-Reactive Protein ELISA Kit (orb219334)

Description	Human C-Reactive Protein ELISA Kit
Reactivity	Human
Range	15.6 pg/ml - 1000 pg/ml
Tested Applications	ELISA
Target	CRP
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Cell Culture Supernatants, Serum, Plasma, Tissue Homogenates Experiment Notes: Colorimetric
Sample Types	Cell Culture Supernates, Serum, Plasma, Tissue Homogenates
Uniprot ID	P02741
Sensitivity	10 pg/ml
Entrez	1401
Expiration Date	Please enquire.



Representative
Standard
Curve