

Product Datasheet

Human Endocan ELISA Kit (orb219357)

Description Sandwich ELISA for quantitative detection of Human Endocan concentrations.

Reactivity Human

Range 15.6 pg/ml - 1000 pg/ml

Tested Applications ELISA

Target ESM1

Storage Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.

Note For research use only

Application notes Cell Culture Supernatants, Serum, Plasma
Experiment Notes: Colorimetric


Sample Types Cell Culture Supernates, Serum, Plasma

Uniprot ID **Q9NQ30**

Sensitivity 7 pg/ml

Entrez **11082**

Expiration Date Please enquire.



Representative
Standard
Curve