


Product Datasheet

Human FAP alpha ELISA Kit (orb219361)

Description	Sandwich ELISA for quantitative detection of Human FAP alpha concentrations.
Reactivity	Human
Range	62.5 pg/ml - 4000 pg/ml
Tested Applications	ELISA
Target	FAP
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Cell Culture Supernatants, Serum, Plasma Experiment Notes: Colorimetric
Sample Types	Cell Culture Supernates, Serum, Plasma
Uniprot ID	Q12884
Sensitivity	10 pg/ml
Entrez	2191
Expiration Date	Please enquire.



Representative
Standard
Curve