

Product Datasheet

Human FR alpha ELISA Kit (orb219376)

Description Sandwich ELISA for quantitative detection of Human FR alpha concentrations.

Reactivity Human

Range 46.9 pg/ml - 3000 pg/ml

Tested Applications ELISA

Target FOLR1

Storage Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.

Note For research use only

Application notes Cell Culture Supernatants, Serum, Plasma, Urine
Experiment Notes: Colorimetric


Sample Types Cell Culture Supernates, Serum, Plasma, Urine

Uniprot ID **P15328**

Sensitivity 10 pg/ml

Entrez **2348**

Expiration Date Please enquire.



Representative
Standard
Curve