


Product Datasheet

Human GH1 ELISA Kit (orb219390)

Description	Sandwich ELISA for quantitative detection of Human GH1 concentrations.
Reactivity	Human
Range	15.6 pg/ml - 1000 pg/ml
Tested Applications	ELISA
Target	GH1
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Cell Culture Supernatants, Serum, Plasma Experiment Notes: Colorimetric
Sample Types	Cell Culture Supernates, Serum, Plasma
Uniprot ID	P01241
Sensitivity	2 pg/ml
Entrez	2688
Expiration Date	Please enquire.



Representative
Standard
Curve