


Product Datasheet

Human GM-CSF ELISA Kit (orb219393)

Description	Sandwich ELISA for quantitative detection of Human GM-CSF concentrations.
Reactivity	Human
Range	15.6 pg/ml - 1000 pg/ml
Tested Applications	ELISA
Target	CSF2
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Cell Culture Supernatants, Serum, Plasma, Urine, Tissue Homogenates Experiment Notes: Colorimetric
Sample Types	Cell Culture Supernates, Serum, Plasma, Tissue Homogenates, Urine
Uniprot ID	P04141
Sensitivity	4 pg/ml
Entrez	1437
Expiration Date	Please enquire.



Representative
Standard
Curve