


Product Datasheet

Human IL12 p40 ELISA Kit (orb219419)

Description	Sandwich ELISA for quantitative detection of Human IL12 p40 concentrations.
Reactivity	Human
Range	62.5 pg/ml - 4000 pg/ml
Tested Applications	ELISA
Target	IL12B
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Cell Culture Supernatants, Serum, Plasma, Tissue Homogenates Experiment Notes: Colorimetric
Sample Types	Cell Culture Supernates, Serum, Plasma, Tissue Homogenates
Uniprot ID	P29460
Sensitivity	30 pg/ml
Entrez	3593
Expiration Date	Please enquire.



Representative
Standard
Curve