


Product Datasheet

Human IL1RA ELISA Kit (orb219433)

Description	Human IL1RA ELISA Kit
Reactivity	Human
Range	31.2 pg/ml - 2000 pg/ml
Tested Applications	ELISA
Target	IL1RN
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Cell Culture Supernatants, Serum, Plasma Experiment Notes: Colorimetric
Sample Types	Cell Culture Supernates, Serum, Plasma
Uniprot ID	P18510
Sensitivity	2 pg/ml
Entrez	3557
Expiration Date	Please enquire.



Representative
Standard
Curve