


Product Datasheet

Human ST2 ELISA Kit (orb219434)

Description	Sandwich ELISA for quantitative detection of human ST2 concentrations.
Reactivity	Human
Range	31.2 pg/ml - 2000 pg/ml
Tested Applications	ELISA
Target	IL1RL1
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Cell Culture Supernatants, Serum, Plasma Experiment Notes: Colorimetric
Sample Types	Cell Culture Supernates, Serum, Plasma
Uniprot ID	Q01638
Sensitivity	15 pg/ml
Entrez	9173
Expiration Date	Please enquire.



Representative
Standard
Curve