


Product Datasheet

Human IL-3 ELISA Kit (orb219445)

Description	Sandwich ELISA for quantitative detection of Human IL-3 concentrations.
Reactivity	Human
Range	31.2 pg/ml - 2000 pg/ml
Tested Applications	ELISA
Target	IL3
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Cell Culture Supernatants, Serum, Plasma Experiment Notes: Colorimetric
Sample Types	Cell Culture Supernates, Serum, Plasma
Uniprot ID	P08700
Sensitivity	15 pg/ml
Entrez	3562
Expiration Date	Please enquire.



Representative
Standard
Curve