


Product Datasheet

Human IL-5 ELISA Kit (orb219451)

Description	Sandwich ELISA for quantitative detection of Human IL-5 concentrations.
Reactivity	Human
Range	1.6 pg/ml - 100 pg/ml
Tested Applications	ELISA
Target	IL5
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Cell Culture Supernatants, Serum, Plasma, Tissue Homogenates Experiment Notes: Colorimetric
Sample Types	Cell Culture Supernates, Serum, Plasma, Tissue Homogenates
Uniprot ID	P05113
Sensitivity	1 pg/ml
Entrez	3567
Expiration Date	Please enquire.



Representative
Standard
Curve