


Product Datasheet

Human IL6 ELISA Kit (orb219452)

Description	Human IL6 ELISA Kit
Reactivity	Human
Range	3.9 pg/ml - 250 pg/ml
Tested Applications	ELISA
Target	IL6
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Cell Culture Supernatants, Serum, Plasma, Tissue Homogenates Experiment Notes: Colorimetric
Sample Types	Cell Culture Supernates, Serum, Plasma, Tissue Homogenates
Uniprot ID	P05231
Sensitivity	2 pg/ml
Entrez	3569
Expiration Date	Please enquire.



Representative
Standard
Curve