


Product Datasheet

Human IL-7 ELISA Kit (orb219454)

Description	Sandwich ELISA for quantitative detection of Human IL-7 concentrations.
Reactivity	Human
Range	7.8 pg/ml - 500 pg/ml
Tested Applications	ELISA
Target	IL7
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Cell Culture Supernatants, Serum, Plasma Experiment Notes: Colorimetric
Sample Types	Cell Culture Supernates, Serum, Plasma
Uniprot ID	P13232
Sensitivity	4 pg/ml
Entrez	3574
Expiration Date	Please enquire.



Representative
Standard
Curve