


Product Datasheet

Human LACI ELISA Kit (orb219465)

Description	Human LACI ELISA Kit
Reactivity	Human
Range	31.2 pg/ml - 2000 pg/ml
Tested Applications	ELISA
Target	TFPI
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Cell Culture Supernatants, Serum, Plasma, UrineExperiment Notes: Colorimetric
Sample Types	Cell Culture Supernates, Serum, Plasma, Urine
Uniprot ID	P10646
Sensitivity	10 pg/ml
Entrez	7035
Expiration Date	Please enquire.



Representative
Standard
Curve