

Product Datasheet

Zika Virus NS1 Antibody (orb2194674)

Description	Zika Virus NS1 Antibody
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length recombinant NS1 protein of Zika virus [Surinarne 71106033] produced in HEK293 cells.
Preservatives	Phosphate Buffered Saline pH 7.4
Form/Appearance	Phosphate Buffered Saline pH 7.4
Storage	Short term: +4C to -2°C. Long Term: -8°C.
Note	For research use only
Application notes	Application Info: Direct ELISA, Western Blot (non-reducing conditions)
Clonality	Monoclonal
Purity	> 95%
Expiration Date	12 months from date of receipt.