


Product Datasheet

Human PDGF-AB ELISA Kit (orb219519)

Description	Human PDGF-AB ELISA Kit
Reactivity	Human
Range	15.6 pg/ml - 1000 pg/ml
Tested Applications	ELISA
Target	PDGFA; PDGFB
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Cell Culture Supernatants, Serum, Plasma Experiment Notes: Colorimetric
Sample Types	Cell Culture Supernates, Serum, Plasma
Uniprot ID	P01127 , P04085
Sensitivity	8 pg/ml
Entrez	5154 , 5155
Expiration Date	Please enquire.



Representative
Standard
Curve