


Product Datasheet

Human Thrombospondin 1 ELISA Kit (orb219552)

Description	Sandwich ELISA for quantitative detection of Human Thrombospondin 1 concentrations.
Reactivity	Human
Range	1560 pg/ml - 100000 pg/ml
Tested Applications	ELISA
Target	THBS1
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Cell Culture Supernatants, Serum, Plasma Experiment Notes: Colorimetric
Sample Types	Cell Culture Supernates, Serum, Plasma
Uniprot ID	P07996
Sensitivity	50 pg/ml
Entrez	7057
Expiration Date	Please enquire.



Representative
Standard
Curve