


Product Datasheet

Human TRK A ELISA Kit (orb219566)

Description	Sandwich ELISA for quantitative detection of Human TRK A concentrations.
Reactivity	Human
Range	156 pg/ml - 10000 pg/ml
Tested Applications	ELISA
Target	NTRK1
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Cell Culture Supernatants, Cell Lysates, Tissue Homogenates Experiment Notes: Colorimetric
Sample Types	Cell Culture Supernates, Cell Lysates, Tissue Homogenates
Uniprot ID	P04629
Sensitivity	10 pg/ml
Entrez	4914
Expiration Date	Please enquire.



Representative
Standard
Curve