

Product Datasheet

Human VEGFR3 ELISA Kit (orb219577)

Description Sandwich ELISA for quantitative detection of Human VEGFR3 concentrations.

Reactivity Human

Range 156 pg/ml - 10000 pg/ml

Tested Applications ELISA

Target FLT4

Storage Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.

Note For research use only

Application notes Cell Culture Supernatants, Serum, Plasma
Experiment Notes: Colorimetric


Sample Types Cell Culture Supernates, Serum, Plasma

Uniprot ID **P35916**

Sensitivity 10 pg/ml

Entrez **2324**

Expiration Date Please enquire.



Representative
Standard
Curve