


Product Datasheet

Mouse VCAM1 ELISA Kit (orb219605)

Description	Mouse VCAM1 ELISA Kit
Reactivity	Mouse
Range	156 pg/ml - 10000 pg/ml
Tested Applications	ELISA
Target	VCAM1
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Cell Culture Supernatants, Serum, Plasma Experiment Notes: Colorimetric
Sample Types	Cell Culture Supernates, Serum, Plasma, Tissue Homogenates
Uniprot ID	P29533
Sensitivity	5 pg/ml
Entrez	22329
Expiration Date	Please enquire.



Representative
Standard
Curve