


Product Datasheet

Mouse CD115 ELISA Kit (orb219606)

Description	Mouse CD115 ELISA Kit
Reactivity	Mouse
Range	62.5 pg/ml - 4000 pg/ml
Tested Applications	ELISA
Target	CSF1R
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Cell Culture Supernatants, Serum, Plasma Experiment Notes: Colorimetric
Sample Types	Cell Culture Supernates, Serum, Plasma
Uniprot ID	P09581
Sensitivity	10 pg/ml
Entrez	12978
Expiration Date	Please enquire.



Representative
Standard
Curve