


Product Datasheet

Mouse CD141 ELISA Kit (orb219613)

Description	Sandwich ELISA for quantitative detection of Mouse CD141 concentrations.
Reactivity	Mouse
Range	62.5 pg/ml - 4000 pg/ml
Tested Applications	ELISA
Target	THBD
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Cell Culture Supernatants, Serum, Plasma, Urine Experiment Notes: Colorimetric
Sample Types	Cell Culture Supernates, Serum, Plasma, Urine
Uniprot ID	P15306
Sensitivity	10 pg/ml
Entrez	21824
Expiration Date	Please enquire.



Representative
Standard
Curve