


Product Datasheet

CTLA4 ELISA kits (orb219615)

Description	CTLA4 ELISA kits
Reactivity	Mouse
Range	156 pg/ml - 10000 pg/ml
Tested Applications	ELISA
Target	CTLA4
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Cell Culture Supernatants, Serum, Plasma Experiment Notes: Colorimetric
Sample Types	Cell Culture Supernates, Serum, Plasma
Uniprot ID	P09793
Sensitivity	10 pg/ml
Entrez	12477
Expiration Date	Please enquire.



Representative
Standard
Curve