


Product Datasheet

Mouse CD268 ELISA Kit (orb219629)

Description	Sandwich ELISA for quantitative detection of Mouse CD268 concentrations.
Reactivity	Mouse
Range	78 pg/ml - 5000 pg/ml
Tested Applications	ELISA
Target	TNFRSF13C
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Cell Culture Supernatants, Serum, Plasma Experiment Notes: Colorimetric
Sample Types	Cell Culture Supernates, Serum, Plasma
Uniprot ID	Q9D8D0
Sensitivity	10 pg/ml
Entrez	72049
Expiration Date	Please enquire.



Representative
Standard
Curve