


Product Datasheet

Mouse CD27 ELISA Kit (orb219631)

Description	Sandwich ELISA for quantitative detection of Mouse CD27 concentrations.
Reactivity	Mouse
Range	62.5 pg/ml - 4000 pg/ml
Tested Applications	ELISA
Target	CD27
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Cell Culture Supernatants, Serum Experiment Notes: Colorimetric
Sample Types	Cell Culture Supernates, Serum
Uniprot ID	P41272
Sensitivity	10 pg/ml
Entrez	21940
Expiration Date	Please enquire.



Representative
Standard
Curve