


Product Datasheet

Mouse CD354 ELISA Kit (orb219637)

Description	Sandwich ELISA for quantitative detection of Mouse CD354 concentrations.
Reactivity	Mouse
Range	31.2 pg/ml - 2000 pg/ml
Tested Applications	ELISA
Target	TREM1
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Cell Culture Supernatants, Serum, Plasma, Tissue Homogenates Experiment Notes: Colorimetric
Sample Types	Cell Culture Supernates, Serum, Plasma, Tissue Homogenates
Uniprot ID	Q9JKE2
Sensitivity	2 pg/ml
Entrez	58217
Expiration Date	Please enquire.



Representative
Standard
Curve