


Product Datasheet

Mouse CD62L ELISA Kit (orb219645)

Description	Sandwich ELISA for quantitative detection of Mouse CD62L concentrations.
Reactivity	Mouse
Range	62.5 pg/ml - 4000 pg/ml
Tested Applications	ELISA
Target	SELL
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Cell Culture Supernatants, Serum, Plasma Notes: Colorimetric
Sample Types	Cell Culture Supernates, Serum, Plasma
Uniprot ID	P18337
Sensitivity	220 pg/ml
Entrez	20343
Expiration Date	Please enquire.



Representative
Standard
Curve