

Product Datasheet

Goat IgG Human Ig Lambda Light Chains (TRITC) antibody (orb21966)

Description
Specificity Reference Reagent

isothiocyanate-conjugated IgG
fraction of polyclo...

Species/Host

Goat

Conjugation

TRITC

Preservatives

Fluorochrome coupled purified goat
IgG lyophilized from a solution in
phosphate buffered saline (pH 7.2).
No preservative added, as it may
interfere with the antibody activity.

Form/Appearance

Fluorochrome coupled purified goat
IgG lyophilized from a solution in
phosphate buffered saline (pH 7.2).
No preservative added, as it may
interfere with the antibody activity.

Storage

Storage: The lyophilized product is
shipped at ambient temperature and
may be stored at +4°C; prolonged
storage at or below -20°C. Prior to
use, an aliquot is thawed slowly in
the dark at ambient temperature,
spun down again and used to
prepare working dilutions by adding
sterile phosphate buffered saline
(PBS, pH 7.2). Repeated thawing and
freezing should be avoided. Working
dilutions should be stored at +4°C,
not refrozen, and preferably used
the same day. If a slight
precipitation occurs upon storage,
this should be removed by
centrifugation. It will not affect the
performance of the product

Note

For research use only

Clonality

Polyclonal

Source

A pool of purified Bence Jones
lambda proteins isolated from
human urine. Freund's complete
adjuvant is used in the first step of
the immunization procedure.

**Hazard
Information**

This product is intended FOR
RESEARCH USE ONLY, and FOR
TESTS IN VITRO, not for use in
diagnostic or therapeutic procedures
involving humans or animals. This

Biorbyt Ltd.

7 Signet Court, Swann's Road, Cambridge, CB5 8LA, United Kingdom

Email: info@biorbyt.com | Phone: +44 (0) 1223 859-353 | Fax: +44 (0)1223 280
240

Biorbyt LLC.

68 TW Alexander Drive
Research Triangle Park
Durham, North
Carolina
27709, United States

Email: info@biorbyt.com | Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558

Expiration Date

12 months from date of receipt.