

Product Datasheet

Mouse E Cadherin ELISA Kit (orb219667)

Description Sandwich ELISA for quantitative detection of Mouse E Cadherin concentrations.

Reactivity Mouse

Range 156 pg/ml - 10000 pg/ml

Tested Applications ELISA

Target CDH1

Storage Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.

Note For research use only

Application notes Cell Culture Supernatants, Serum, Plasma
Experiment Notes: Colorimetric


Sample Types Cell Culture Supernates, Serum, Plasma

Uniprot ID **P09803**

Sensitivity 15 pg/ml

Entrez **12550**

Expiration Date Please enquire.



Representative
Standard
Curve