


## Product Datasheet

### Mouse Galectin 3 ELISA Kit (orb219679)

|                            |                                                                                               |
|----------------------------|-----------------------------------------------------------------------------------------------|
| <b>Description</b>         | Mouse Galectin 3 ELISA Kit                                                                    |
| <b>Reactivity</b>          | Mouse                                                                                         |
| <b>Range</b>               | 62.5 pg/ml - 4000 pg/ml                                                                       |
| <b>Tested Applications</b> | ELISA                                                                                         |
| <b>Target</b>              | LGALS3                                                                                        |
| <b>Storage</b>             | Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles. |
| <b>Note</b>                | For research use only                                                                         |
| <b>Application notes</b>   | Cell Culture Supernatants, Serum, Plasma<br>Experiment Notes: Colorimetric                    |
| <b>Sample Types</b>        | Cell Culture Supernates, Serum, Plasma                                                        |
| <b>Uniprot ID</b>          | <b>P16110</b>                                                                                 |
| <b>Sensitivity</b>         | 10 pg/ml                                                                                      |
| <b>Entrez</b>              | <b>16854</b>                                                                                  |
| <b>Expiration Date</b>     | Please enquire.                                                                               |



Representative  
Standard  
Curve