


Product Datasheet

Mouse HVEM ELISA Kit (orb219687)

Description	Sandwich ELISA for quantitative detection of Mouse HVEM concentrations.
Reactivity	Mouse
Range	15.6 pg/ml - 1000 pg/ml
Tested Applications	ELISA
Target	TNFRSF14
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Cell Culture Supernatants, Serum, Plasma Experiment Notes: Colorimetric
Sample Types	Cell Culture Supernates, Serum, Plasma
Uniprot ID	Q80WM9
Sensitivity	1 pg/ml
Entrez	230979
Expiration Date	Please enquire.



Representative
Standard
Curve