


Product Datasheet

Mouse IL4 ELISA Kit (orb219724)

Description	Mouse IL4 ELISA Kit
Reactivity	Mouse
Range	31.2 pg/ml - 2000 pg/ml
Tested Applications	ELISA
Target	IL4
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Cell Culture Supernatants, Serum, Plasma, Tissue Homogenates Experiment Notes: Colorimetric
Sample Types	Cell Culture Supernates, Serum, Plasma, Tissue Homogenates
Uniprot ID	P07750
Sensitivity	15 pg/ml
Entrez	16189
Expiration Date	Please enquire.



Representative
Standard
Curve