


Product Datasheet

Mouse M-CSF ELISA Kit (orb219740)

Description	Sandwich ELISA for quantitative detection of Mouse M-CSF concentrations.
Reactivity	Mouse
Range	15.6 pg/ml - 1000 pg/ml
Tested Applications	ELISA
Target	CSF1
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Cell Culture Supernatants, Serum, Plasma, Tissue Homogenates Experiment Notes: Colorimetric
Sample Types	Cell Culture Supernates, Serum, Plasma, Tissue Homogenates
Uniprot ID	P07141
Sensitivity	1 pg/ml
Entrez	12977
Expiration Date	Please enquire.



Representative
Standard
Curve