


Product Datasheet

Mouse MERTK ELISA Kit (orb219741)

Description	Mouse MERTK ELISA Kit
Reactivity	Mouse
Range	78 pg/ml - 5000 pg/ml
Tested Applications	ELISA
Target	MERTK
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Cell Culture Supernatants, Cell Lysates, Tissue Homogenates, SerumExperiment Notes: Colorimetric
Sample Types	Cell Culture Supernates, Serum, Cell Lysates, Tissue Homogenates
Uniprot ID	Q60805
Sensitivity	10 pg/ml
Entrez	17289
Expiration Date	Please enquire.



Representative
Standard
Curve