


Product Datasheet

Mouse NGAL ELISA Kit (orb219750)

Description	Mouse NGAL ELISA Kit
Reactivity	Mouse
Range	15.6 pg/ml - 1000 pg/ml
Tested Applications	ELISA
Target	LCN2
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Cell Culture Supernatants, Serum, Plasma, Urine Experiment Notes: Colorimetric
Sample Types	Cell Culture Supernates, Serum, Plasma, Urine
Uniprot ID	P11672
Sensitivity	8 pg/ml
Entrez	16819
Expiration Date	Please enquire.



Representative
Standard
Curve