


Product Datasheet

Mouse LOX-1 ELISA Kit (orb219752)

Description	Mouse LOX-1 ELISA Kit
Reactivity	Mouse
Range	31.2 pg/ml - 2000 pg/ml
Tested Applications	ELISA
Target	OLR1
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Cell Culture Supernatants, Serum, Plasma Experiment Notes: Colorimetric
Sample Types	Cell Culture Supernates, Serum, Plasma
Uniprot ID	Q9EQ09
Sensitivity	1 pg/ml
Entrez	108078
Expiration Date	Please enquire.



Representative
Standard
Curve