


Product Datasheet

Mouse Thrombopoietin ELISA Kit (orb219776)

Description	Mouse Thrombopoietin ELISA Kit
Reactivity	Mouse
Range	31.2 pg/ml - 2000 pg/ml
Tested Applications	ELISA
Target	THPO
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Cell Culture Supernatants, Serum, Plasma Experiment Notes: Colorimetric
Sample Types	Cell Culture Supernates, Serum, Plasma
Uniprot ID	P40226
Sensitivity	5 pg/ml
Entrez	21832
Expiration Date	Please enquire.



Representative
Standard
Curve