


Product Datasheet

Mouse TIE2 ELISA Kit (orb219777)

Description	Mouse TIE2 ELISA Kit
Reactivity	Mouse
Range	125 pg/ml - 8000 pg/ml
Tested Applications	ELISA
Target	TEK
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Cell Culture Supernatants, Serum, Plasma, Tissue Homogenates Experiment Notes: Colorimetric
Sample Types	Cell Culture Supernates, Serum, Plasma, Tissue Homogenates
Uniprot ID	Q02858
Sensitivity	10 pg/ml
Entrez	21687
Expiration Date	Please enquire.



Representative
Standard
Curve