


Product Datasheet

Mouse YKL-40 ELISA Kit (orb219788)

Description	Sandwich ELISA for quantitative detection of Mouse YKL-40 concentrations.
Reactivity	Mouse
Range	62.5 pg/ml - 4000 pg/ml
Tested Applications	ELISA
Target	CHI3L1
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Cell Culture Supernatants, Serum, Plasma Experiment Notes: Colorimetric
Sample Types	Cell Culture Supernates, Serum, Plasma
Uniprot ID	Q61362
Sensitivity	10 pg/ml
Entrez	12654
Expiration Date	Please enquire.



Representative
Standard
Curve