


Product Datasheet

Rat CD62L ELISA Kit (orb219799)

Description	Sandwich ELISA for quantitative detection of Rat CD62L concentrations.
Reactivity	Rat
Range	625 pg/ml - 40000 pg/ml
Tested Applications	ELISA
Target	SELL
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Cell Culture Supernatants, Serum, Plasma Experiment Notes: Colorimetric
Sample Types	Cell Culture Supernates, Serum, Plasma
Uniprot ID	P30836
Sensitivity	20 pg/ml
Entrez	29259
Expiration Date	Please enquire.



Representative
Standard
Curve