


Product Datasheet

Rat IL12 p70 ELISA Kit (orb219824)

Description	Sandwich ELISA for quantitative detection of Rat IL-12 p70 concentrations.
Reactivity	Rat
Range	15.6 pg/ml - 1000 pg/ml
Tested Applications	ELISA
Target	IL12A; IL12B
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Cell Culture Supernatants, Serum, Plasma, Tissue Homogenates Experiment Notes: Colorimetric
Sample Types	Cell Culture Supernates, Serum, Plasma, Tissue Homogenates
Uniprot ID	Q9R278, Q9R103
Sensitivity	15 pg/ml
Entrez	64546, 84405
Expiration Date	Please enquire.



Representative
Standard
Curve