


Product Datasheet

Rat IL-2 ELISA Kit (orb219828)

Description	Sandwich ELISA for quantitative detection of Rat IL-2 concentrations.
Reactivity	Rat
Range	31.2 pg/ml - 2000 pg/ml
Tested Applications	ELISA
Target	IL2
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Cell Culture Supernatants, Serum, Plasma, Tissue Homogenates Experiment Notes: Colorimetric
Sample Types	Cell Culture Supernates, Serum, Plasma, Tissue Homogenates
Uniprot ID	P17108
Sensitivity	7 pg/ml
Entrez	116562
Expiration Date	Please enquire.



Representative
Standard
Curve