


Product Datasheet

Rat IL6 ELISA Kit (orb219833)

Description	Rat IL6 ELISA Kit
Reactivity	Rat
Range	62.5 pg/ml - 4000 pg/ml
Tested Applications	ELISA
Target	IL6
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Cell Culture Supernatants, Serum, Plasma, Tissue Homogenates Experiment Notes: Colorimetric
Sample Types	Cell Culture Supernates, Serum, Plasma, Tissue Homogenates
Uniprot ID	P20607
Sensitivity	30 pg/ml
Entrez	24498
Expiration Date	Please enquire.



Representative
Standard
Curve