


Product Datasheet

Rat Inhibin beta A ELISA Kit (orb219834)

Description	Sandwich ELISA for quantitative detection of Rat Inhibin beta A concentrations.
Reactivity	Rat
Range	15.6 pg/ml - 1000 pg/ml
Tested Applications	ELISA
Target	INHBA
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Cell Culture Supernatants, Serum, Plasma Experiment Notes: Colorimetric
Sample Types	Cell Culture Supernates, Serum, Plasma
Uniprot ID	P18331
Sensitivity	12 pg/ml
Entrez	29200
Expiration Date	Please enquire.



Representative
Standard
Curve