


Product Datasheet

Rat MMP2 ELISA Kit (orb219840)

Description	Sandwich ELISA Kit of Rat MMP2
Reactivity	Rat
Range	625 pg/ml - 40000 pg/ml
Tested Applications	ELISA
Target	MMP2
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Cell Culture Supernatants, Serum, Plasma Experiment Notes: Colorimetric
Sample Types	Cell Culture Supernates, Serum, Plasma
Uniprot ID	P33436
Sensitivity	300 pg/ml
Entrez	81686
Expiration Date	Please enquire.



Representative
Standard
Curve