


Product Datasheet

Rat NGAL ELISA Kit (orb219842)

Description	Rat NGAL ELISA Kit
Reactivity	Rat
Range	78 pg/ml - 5000 pg/ml
Tested Applications	ELISA
Target	LCN2
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Cell Culture Supernatants, Serum, Plasma, UrineExperiment Notes: Colorimetric
Sample Types	Cell Culture Supernates, Serum, Plasma, Urine
Uniprot ID	P30152
Sensitivity	10 pg/ml
Entrez	170496
Expiration Date	Please enquire.



Representative
Standard
Curve