


Product Datasheet

Rat PDGF-AB ELISA Kit (orb219844)

Description	Rat PDGF-AB ELISA Kit
Reactivity	Rat
Range	31.2 pg/ml - 2000 pg/ml
Tested Applications	ELISA
Target	PDGFA; PDGFB
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Cell Culture Supernatants, Serum, Plasma Experiment Notes: Colorimetric
Sample Types	Cell Culture Supernates, Serum, Plasma
Uniprot ID	P28576, Q05028
Sensitivity	3 pg/ml
Entrez	25266, 24628
Expiration Date	Please enquire.



Representative
Standard
Curve