



Product Datasheet

Acid Phosphatase Assay Kit (orb219857)

Description	Acid Phosphatase Assay Kit
Reactivity	All
Range	0.05 mmol/L - 5 mmol/L
Tested Applications	FA
Storage	Shipped and store at 4°C for 6 months.
Note	For research use only
Application notes	Experiment Notes: Colorimetric
Sample Types	Urine, Serum, Plasma, Tissue extracts, Cell lysate, Cell culture media, Other biological fluids
CAS Number	9001-77-8
Expiration Date	12 months from date of receipt.



Representative Standard Curve



Positive Control reaction in 96-well pla...