


Product Datasheet

Alanine Transaminase Assay Kit (orb219858)

Description	Alanine Transaminase Assay Kit
Reactivity	All
Range	0.625 µmol/ml - 20 µmol/ml
Tested Applications	FA
Storage	Shipped and store at 4°C for 6 months.
Note	For research use only
Application notes	Experiment Notes: Colorimetric
Sample Types	Urine, Serum, Plasma, Tissue extracts, Cell lysate, Cell culture media, Other biological fluids
CAS Number	9000-86-6
Expiration Date	12 months from date of receipt.



Representative
Standard
Curve