


Product Datasheet

Glutamate Assay Kit (orb219861)

Description	Glutamate Assay Kit
Reactivity	All
Range	0.03 mmol/L - 3 mmol/L
Tested Applications	FA
Storage	Shipped and store at 4°C for 6 months.
Note	For research use only
Application notes	Experiment Notes: Colorimetric
Sample Types	Urine, Serum, Plasma, Tissue extracts, Cell lysate, Cell culture media, Other biological fluids
CAS Number	56-86-0
Expiration Date	12 months from date of receipt.



Representative
Standard
Curve